

ReverTra Ace

(MMLV Reverse Transcriptase RNaseH-)

Code No. TRT-

Lot No. *****

Storage Store at -20°C

Size 2,000units(101T) 10,000units(101)

Source : *Escherichia coli* DH5 α (pDV68)

Concentration : 100 units/ μ l

Unit Definition : One unit is the amount of enzyme activity that incorporates 1 nmole of TTP into acid insoluble fraction in 10 minutes at 42°C when poly(A)⁺ RNA and oligo (dT)₂₀ are used as template-primer.

Assay Condition : 50 mM Tris-HCl(pH8.3)
75 mM KCl
3 mM MgCl₂
10 mM DTT
240 μ g/ml poly(A)⁺ RNA
0.1 μ M oligo (dT)₂₀
1 mM dTTP
0.3 μ M [methyl-³H]TTP

Storage Buffer : 20 mM Tris-HCl(pH7.5)
100 mM NaCl
0.1 mM EDTA
1 mM DTT
0.01 % Nonidet P-40
50 % Glycerol

5 X Buffer : Not containing dNTPs. cDNA synthesis reaction requires 1mM of dNTPs.

Quality Control

1. RNase Activity : When 50 units of this enzyme were incubated with 0.5 μ g of MS2 RNA for 1 hour at 42°C, RNase activity was not observed after agarose gel electrophoresis.

Remarks Research Use Only. Not for Human Use.