Code No. FIK-203  
保存温度 -20℃

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## Label License Statement

### **NOTICE TO PURCHASER: LIMITED LICENSE**

This product is provided under an Agreement between Biotium and TOYOBO. The manufacture, use, sale, offer for sale, or import of this product is subject to one or more patents or pending applications owned or licensed by Biotium, Inc. (U.S. Patent Nos. US7,803,943 B2, US7,776, 567 B2, Japan Patent No. 5249013). The purchase price of this product includes a limited, non-transferable immunity from suit under U.S. Patent Nos. US7,803,943 B2, US7,776, 567 B2, Japan Patent No. 5249013 and corresponding patent claims outside of the United States, to use only the purchased amount of the product and components of the product solely for the buyer's own internal research (whether the buyer is an academic or for-profit entity), in a manner consistent with the accompanying product literature. No right to use, sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party, or otherwise use this product or its components or materials made using this product or its components for any other purpose or a Commercial Purpose, is hereby granted expressly, by implication or by estoppel. Commercial Purpose means any activity by a party for consideration and include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service including without limitation, reporting information or data resulting from use of the product or a its components for a fee or other commercial consideration; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. For products that are subject to multiple limited use label licenses, the most restrictive terms apply. This product is for research use only. Diagnostic uses require a separate license from Biotium, Inc. For information on purchasing a license to this product including for purposes other than research, contact Director of Licensing, Biotium, Inc., 3159 Corporate Place, Hayward, CA 94545, Tel: +1-(510) 265-1027. Fax: +1-(510) 265-1352. Email: [btinfo@biotium.com](mailto:btinfo@biotium.com). All names containing the designation® are trademarks of Biotium.

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