

Phage Lambda DNA

Code No. DNA-001
Lot No. *****
Storage Store at -20
Size 400 µg

C e r t i f i c a t e o f A n a l y s i s

Product : Deoxyribonucleic acid -DNA

This vector is modified by methylase, because it was produced by E.coli W3350. Some restriction enzyme that are affected by methylation can not digest it.

Source : Bacteriophage c 857 Sam 7

Appearance : Solution in 10mM Tris-HCl pH8.0, 1mM EDTA

Concentration : ***** O.D.DNA / ml

Volume : ***** ml / tube

O.D.260 / O.D.280 : *****

Molecular Weight / Size : 32 Mdal , 48502 bps