

KOD One™ PCR Master Mix -Blue-Code No. **KMM-201**Lot No. *********Size **1mL × 5**

Content : **KOD One™ PCR Master Mix -Blue-** **1mL × 5**
(200 reactions,
50 μ l total reaction volume)

Appearance : **Blue aqueous solution**

Quality Control

Functional test

The 10 kb fragment could be amplified from human genomic DNA using fast cycle which extension time is 50 sec.



manual

TOYOBO CO.,LTD.

Biotechnology Operating Department

Osaka Umeda Twin Towers South,
1-13-1 Umeda, Kita-ku, Osaka 530-0001, Japan
Tel:+81-6-6348-3846(Global), 3786(Japan)
Fax:+81-6-6348-3833
www.toyobo-global.com/products/bio/