

ΦX174/HaeIII digest

Code No. DNA-012
Lot No. *****
Storage Store at -20°C
Size 15 µg

Certificate of Analysis

Product : ΦX174 / *Hae*III digest
A *Hae*III digest of ΦX174 am 3 cs70 RFI DNA was dissolved in 10mM Tris-HCl and 1mM EDTA.

Concentration : * * * * µg/µl

Volume : * * * * µl/tube(15µg/tube)

Fragment sizes :

fragment No.	Number of Base Pairs(bps)	Molecular weight ($\times 10^5$ daltons)
1	1,353	8.93
2	1,078	7.11
3	872	5.76
4	603	3.98
5	310	2.05
6B	281	1.85
6A	271	1.79
7	234	1.54
8	194	1.28
9	118	0.78
10	72	0.48

Handling recommendations : 1. Centrifuge tube before opening to improve recovery of content.