

## ΦX174/HaeIII digest

Code No. DNA-012

Lot No. \*\*\*\*\*

Size 15 μg

### Certificate of Analysis

Product : ΦX174 / HaeIII digest  
 A HaeIII digest of ΦX174 am 3 cs70 RFI DNA was dissolved in 10mM Tris-HCl and 1mM EDTA.

Concentration :       \* \* \* \*       μg/μLVolume :       \* \* \* \*       μL/tube(15μg/tube)

Fragment sizes	fragment No.	Number of Base Pairs(bps)	Molecular weight (× 10 <sup>5</sup> daltons)
	1	1,353	8.93
	2	1,078	7.11
	3	872	5.76
	4	603	3.98
	5	310	2.05
	6B	281	1.85
	6A	271	1.79
	7	234	1.54
	8	194	1.28
	9	118	0.78
	10	72	0.48

Handling recommendations : 1. Centrifuge tube before opening to improve recovery of content.

