

Loading Quick™

ΦX174/HaeIII digest

Code No. DNA-112

Lot No. *****

Size 20 µg

Product : A HaeIII digest of ΦX174 am 3 cs70 RFI DNA was dissolved in 10mM Tris-HCl, 15mM NaCl, 8mM EDTA, 0.0125% Brom Phenol Blue and 3.75% Ficoll.

Concentration : 0.04 µg/µL

Volume : 500 µL/tube(20µg/tube)

Fragment sizes	fragment No.	Number of Base Pairs(bps)	Molecular weight (× 10 ⁵ daltons)
	1	1,353	8.93
	2	1,078	7.11
	3	872	5.76
	4	603	3.98
	5	310	2.05
	6B	281	1.85
	6A	271	1.79
	7	234	1.54
	8	194	1.28
	9	118	0.78
	10	72	0.48

Handling recommendations : Centrifuge tube before opening to improve recovery of content.

