

KOD One™ PCR Master Mix

Code No. **KMM-101**

Lot No. *********

Size **1 mL x 5**

Content : **KOD One™ PCR Master Mix** 1 mL x 5
(200 reactions,
50 µl total reaction volume)

Appearance : **Colorless and transparent aqueous solution**

Quality Control

Functional test

The 10 kb fragment could be amplified from human genomic DNA using fast cycle which extension time is 50 sec.



manual

TOYOBO CO., LTD.

Life Science Department

2-8, Dojima Hama 2-chome, Kita-ku, Osaka, 530-8230, Japan
TEL : 06-6348-3786 FAX : 06-6348-3833
Technical Service : 06-6348-3888
【URL】 <http://lifescience.toyo.co.jp/>