

KOD One™ PCR Master Mix -Blue-

Code No. KMM-201

Lot No. *****

Size 1 mL x 5

Content : KOD One™ PCR Master Mix -Blue- 1 mL x 5
(200 reactions,
50 μ l total reaction volume)

Appearance : Blue aqueous solution

Quality Control

Functional test

The 10 kb fragment could be amplified from human genomic DNA using fast cycle which extension time is 50 sec.



manual

TOYOBO CO., LTD.

Life Science Department

2-8, Dojima Hama 2-chome, Kita-ku, Osaka, 530-8230, Japan
TEL : 06-6348-3786 FAX : 06-6348-3833
Technical Service : 06-6348-3888
【URL】 <http://lifescience.toyo.co.jp/>