

RNase Inhibitor, Recombinant

Code No. SIN-201

Lot No. *****

Size 2,500units

Source : *E.coli* cells expressing a recombinant clone.

Concentration : 40 units/ μ L

Unit Definition : One unit is defined as the amount of RNase inhibitor, Recombinant required to inhibit the activity of 5 ng of ribonuclease A by 50 %. Activity is measured by the inhibition of hydrolysis of cytidine 2', 3'-cyclic monophosphate by ribonuclease A

Storage Buffer : 20 mM HEPES-KOH (pH7.6)
50 mM KCl
8 mM DTT
50 % Glycerol
<0.01 % Tween™ 20

Usage Notes : Since ribonucleases typically retain activity under denaturing conditions, care must be taken to avoid denaturing RNase inhibitor molecules which have complexed with ribonuclease. To prevent the release of active ribonuclease, temperature greater than 50°C and high concentration of urea or other denaturing agents should be avoided.
RNase inhibitor, Recombinant is active over a broad pH range.
Use ~1unit of inhibitor per μ L of solution.

References

1. Blackburn, P. and Moore, S. (1982) In: *The Enzymes* Vol. XV, Part B, Academic Press, NY.
2. Blackburn, P., Wilson, G. and Moore, S. (1977) Ribonuclease inhibitor from human placenta. Purification and Properties. *J. Biol. Chem.* **252**, 5904.

All trademarks, trade names, or company names referenced herein are used for identification only and are the property of their respective owners.



information

TOYOBO CO.,LTD.

Biotechnology Operating Department

Osaka Umeda Twin Towers South,
1-13-1 Umeda, Kita-ku, Osaka 530-0001, Japan
Tel:+81-6-6348-3846(Global), 3786(Japan)
Fax:+81-6-6348-3833
www.toyobo-global.com/products/bio/